



## Titrator Compact G10S

### **Robust titrator for general applications.**

*The G10S introduces the solid, simple operation and space saving design of the Compact titrators. Up to 5 analysis methods with 10 method functions + 4 shortcuts per user. Ideal for dedicated QC applications with several operators.*



### **Predefined Methods**

Includes 20 optimized methods for common applications in all industries. Use directly or easily edit them to suit your needs.

### **One Click User Interface**

With the intuitive One Click touchscreen interface, customize home screens to run your workflows quickly and easily with a single tap on the screen.

### **Plug & Play Accessories**

Sensors and burettes are automatically recognized when connected. Avoid errors by eliminating extra data entry steps and reduce instrument setup time.

*Material No.: 30267117*

[Overview](#)[Documentation](#)[Accessories](#)[Consumables](#)

## Specifications - Titrator Compact G10S

Max. number of methods	5
Method Functions per Method	10 max.
Additional Burette Drives (Optional)	None
LabX Software	Not compatible
Shortcuts per user	4
Application Type	Potentiometric
Measurement Type	Chromatometric Amperometric indication Diazotitration Redox Complexometry Direct measurement (ion selective) Potentiometric Precipitation Voltametric indication Photometric Cerimetric Bromatometric Acid/base Arsenometric
Karl Fischer titration	No
User Rights and Management	Predefined Expert" and "Routine" roles"
Solvent Manager Option	Yes
Parallel Titration	No
Samples per Analysis (Max.)	120
Barcode Reader Connection	No
Predefined methods	20


[Online Chat](#)

User-Defined Calculations	No
Task List (Max.)	0
One Click® User Interface	Yes
Flexibility and Modularity	Automation Plug-and-play accessories Sample weight transfer SmartChemicals PC software One-click user interface User-defined calculations Pre-installed methods User management
Plug & Play Features	Burettes and Sensors
Compliance and Data Security	User management (fixed groups)
Printer Connection	Yes (USB)
Connectivity and Software	Network export USB export
Automation Option(s)	Rondolino
Display	Full-color 5.7 in VGA, touch screen
Languages	Chinese, English, French, German, Italian, Korean, Polish, Portuguese, Russian, Spanish



Online Chat

## Secure Operation

Designate "Expert" and "Routine" user groups to allow full or limited access to methods, results and instrument setup parameters.



## Simple Automation

Compact titrators support an automatic Rondolino sample changer. Analysis of up to 9 samples plus automated cleaning and conditioning functions speed up your daily tasks.



## Safe Chemical Handling

Reduce user contact with chemicals by automatically draining and refilling the titrator cell and preventing waste overflow with dedicated accessories.



## Flexible Data Management

Compact titrators offer multiple options for data handling. Export data in CSV, PDF or XML files via a memory stick or connect a printer directly via the USB port.



## Documentation

### Product Brochures



## Compact Titrators

Taking the step from manual to fully automated titration allows you to increase efficiency by dedicating your time to other important tasks. With its...

---

▼ Datasheets

---

▼ UserCom

---

▼ Applications

---

^ Accessories



### Bottle Adapter (Merck)

Merck bottle adapter

*Material No.: 23774*

**USD 74.70**

Add to Cart



Online Chat

### Brown Glass Bottle, 1L

Brown glass bottle, 1L

Material No.: 71296

Get a Quote



### Cable RS9-RS9 (M/W) 1m

RS232 9-pin male to female cable for connecting a balance to a printer, PC or titrator

Material No.: 11101051

**USD 145.00**

Add to Cart



Online Chat

### Cable USB (1.8 m A-B) for PC or printer

USB cable, 1.8 m

Material No.: 51191926

**USD 35.50**

Add to Cart



### Cable USB 412

Material No.: 21901309

Get a Quote



## DV1000 interchangeable burettes

The DV1000 series of interchangeable burettes allow secure and rapid titrant change. With the choice of four different burette sizes you can optimize the dispensing volume and achieve highest precision.

Material No.: 51107503

Get a Quote



## DV704 Titration Stand



The external **DV704** titration stand can be used for volumetric and coulometric Karl Fischer titrations and for potentiometric titrations. The titration stand includes a magnetic stirrer and can be employed with all Excellence and Compact titrators from METTLER TOLEDO.

The **DV704** can be optimally combined with the Solvent Manager to change solvent easily, safely, and without operator contact, for example for the Karl Fischer titration. All the usual titration accessories such as sensors, standard beakers or overhead stirrer can of course be directly utilized with the **DV704** or with the external titration arm.

Material No.: 51109259

Get a Quote



## External power supply 100 - 240 VAC

Power supply

Material No.: 51105795

**USD 201.00**

Add to Cart



## External SmartChemical Reader

Material No.: 30486743

Get a Quote



Online Chat

## Fisher Bottle Adapter

Fisher bottle adapter

Material No.: 23787

**USD 64.00**

Add to Cart





## Magnetic Stirrer Bar

Magnetic stirrer bar

Material No.: 51191159

**USD 25.90**

Add to Cart



## Sample Changer Rondolino

Automated Titration Stand

Material No.: 51108500

Get a Quote



Online Chat

## Solvent Manager

The Solvent Manager is a device for the handling of liquids against the contact of the user with chemicals.

Material No.: 51105652

Get a Quote

[Consumables](#)

[Accessories](#)



## Beaker PP (100mL), 480 pcs.

Disposable titr. beakers 480

Material No.: 51109389

**USD 298.00**

[Add to Cart](#)

- ✓ EasyPlus™ Accessories
- ✓ Metal sensors
- ✓ Photometric sensor
- ✓ Surfactant sensors
- ✓ Temperature probe and reference sensors
- ✓ pH glass sensors



Online Chat

**Get Your Quote**  
**Titration Compact G10S**

[Request Your Quote Today](#)



## Your Recently Viewed Products and Services

[EasyPlus Titrator Easy EPM](#)

[EasyPlus Dosing Device Easy D...](#)

[Sodium Analyzer Easy Na](#)

---

### Products & Solutions

[Laboratory Balances](#)

[Industrial Scales](#)

[Retail Weighing Scales](#)

[Rainin Pipettes and Tips](#)

[Analytical Instruments](#)

[Automated Reactors and In Situ Analysis](#)

[Process Analytics](#)

[Inline Product Inspection](#)

[Transport and Logistics](#)



Online Chat

---

### Services & Support

[Service Offerings](#)

[Product Registration](#)

[Technical Support](#)

[Payment Methods](#)

---

### Events & Expertise

## [Expertise Library](#)

[Literature: White Papers, Guides, Brochures](#)

[Technical Documentation](#)

[On Demand Webinars](#)

[Videos](#)

## [Live Events](#)

[Seminars & Training](#)

[Trade Shows & Exhibitions](#)

[Live Webinars](#)

## **About Us**

---

### [Jobs & Careers](#)

METTLER TOLEDO

[Sustainability](#)

[Business Description](#)

[Digitalization](#)

[Global Presence](#)

[Innovation & Quality](#)

[Management](#)

[Investor Relations](#)

[Newsroom](#)



Online Chat

## **Quick Links**

---

[Contact Us](#)

[Search](#)

[Document Finder](#)

[Product Finder](#)

[Service Finder](#)

[Videos](#)

---

© METTLER TOLEDO. All rights reserved.

[Do not sell my personal information](#)

[Legal, Terms, Privacy](#)

[Contact us](#)

[Site Map](#)

[Global Presence](#)

---



Online Chat